

Wash-Checks H & Lumen Devices

Valisafe[®]

clearly a better indicator



Easy to
use

- For monitoring of the cleaning efficacy when cleaning rigid and flexible lumen devices.
- Test device is designed to mimic the removal of blood and tissue from inside lumen devices
- Test soil contains no blood components and is non toxic
- Easy to read
- V3503935 Indicator fits universally into both the H and E wash check holders
- 2 year shelf life
- Durable stainless steel holder

Valisafe Wash-Checks H & Lumen Devices

The Valisafe wash-check cleaning efficacy system for rigid and flexible lumen instruments is designed to monitor the cleaning efficacy of your flow and pressure systems used to clean lumens during the wash phase in your washer disinfectant. Used at least twice a day it provides routine evidence of the efficacy of your processes. Successfully clean monitors can then be stored and maintained as part of your record system.

Order Details

- V3503935 Wash-Check cleaning indicators (Wash check H)
25 indicators per vial - 4 vials/box
- V3503940 Wash-Check cleaning indicator holder – Rigid lumen (Wash-Check H holder)
1 holder
- V3503945 Wash-Check cleaning indicator holder – flexible lumen (Wash-Check E holder)
1 holder



V3503935 Wash-Check cleaning indicators (Wash check H)



V3503940 Wash-Check cleaning indicator holder - Rigid lumen



V3503945 Wash-Check cleaning indicator holder

Name:	Wash-Check H lumen cleaning indicator Wash-Check H indicator holder - rigid lumen device Wash-Check E indicator holder flexible lumen device
Classification:	AAMI ISO 15883
Storage conditions:	room temperature in sealed foil bag
Lot no:	yes
Prod date:	yes
Exp date:	yes
Correct symbols:	In accordance with AAMI ISO 15883
Indicator:	Cleaning efficacy indicator
Non toxic:	yes
Disp/single use::	indicator single use - holder multiple uses