

Class 4 Steam Indicator

Valisafe®

clearly a better indicator

- Conforms to EN ISO 11140-1 Class 4 / AAMI ST79
- Non toxic
- Easy to use
- Recordable
- Clear colour change pink to brown
- 3 year shelf life



Class 4 Steam Indicator

The Valisafe Class 4 indicator is designed to function in a steam sterilizer operating at 21°C(250°F) -134°C (273°F) for 3.5 minutes or longer. The steam indicator strip gives clear visible confirmation that the critical parameters have been attained at the location of the indicator strip within the sterilizer chamber.

Order Details

V3503975

Class 4 Steam Indicator

Bags of 25 indicators, 24 bags in a case

STEAM INDICATOR

STEAM
LOT
000000







ISO 11140-1 Class 4
 Pink/Rosa/Rosa = Underexposed/Insuficiente/Insuficiente
 Brown/Marrone/Marrón = Correct/Correcto/Correcto



Art No.: V3503975

Name:	Class 4 Steam Indicator
Classification:	EN ISO 11140-1 Class 4 and AAMI ST79
Operating temp:	21°C (250°F) -134°C (273°F)
Time:	3.5 minutes
Storage conditions:	50-100°F (10-38°C)
RH:	10%-90%
Lot no.	yes
Prod date:	on the packaging
Exp date:	3 years
Correct symbols:	In accordance with EN ISO 11140-1
Indicator:	changes from pink to dark brown
Disp/single use:	yes
Physical Charateristics:	
Dimensions	0.5625" w x 3.8"L
Coated	yes
Extra information:	
Non toxic/leadfree	yes

Before
Original indicator



After
After total process

