

Dual Steam Indicator

Valisafe®

clearly a better indicator

- Conforms to EN ISO 11140-1 Class 4 / AAMI ST79
- Non toxic
- Easy to use
- Recordable
- Clear colour change pink to brown / black
- 3 year shelf life



Dual Steam Indicator

The Valisafe Dual Steam Indicator is a multi-variable indicator. It gives visual confirmation of the attainment of the critical parameters inside a steam sterilizer operating at 121°C(250°F) -134°C (273°F) for 3.5 minutes or longer at the location of the indicator within the steriliser chamber.

Order Details

V3503967

Dual Steam Indicator

Bags of 250, 16 bags in a case

Name:	Dual Steam Indicator
Classification:	EN ISO 11140-1 Class 4 and AAMI ST79
Operating temp:	121°C (250°F) – 134°C (273°)
Time:	3.5 minutes
Storage conditions:	50-100°F (10-38°C)
RH:	10%-70%
Lot no.	yes
Prod date:	on box
Exp date:	3 years
Correct symbols:	In accordance with EN ISO 11140-1 Class 4
Indicator:	changes from pink to dark brown
Disp/single use:	yes
Physical Characteristics:	
Dimensions	0.5635 "w x 7.66 " L
Usage	The strip can be cut in 2 and each end used as an individual strip
Extra information:	
Non toxic/leadfree	yes
Recordable.	yes

DUAL MULTI-PARAMETER INDICATOR

Class 4
121 °C - 134 °C
Art. No.: V3502038

Dual multi-parameter indicator
EN-ISO 11140-1
Sterilizing conditions
turn indicator from
pink to brown/black.
LOT 100305

Valisafe
STEAM

DUAL MULTI-PARAMETER INDICATOR

Class 4
121 °C - 134 °C
Art. No.: V3502038

Dual multi-parameter indicator
EN-ISO 11140-1
Sterilizing conditions
turn indicator from
pink to brown/black.
LOT 100305

Valisafe
STEAM



Before
Original indicator



After
After total process

